

**Supplemental Table 2.** Accuracy studies for Roche performed using a recovery test for 25(OH)D (nmol/L)

	Background value (nmol/L, <i>n</i> = 2)	Additive concentration (nmol/L)	Measured concentration (nmol/L, <i>n</i> = 2)	Recovery rate (%)
Centaur	36.7	52.3	37	75.9
		315	49.3	51.7
Elecsys	41.1	52.3	41.3	82.4
		315	56.8	62.9
LC-MS/MS	62.9	52.3	63.1	124.1
		315	86.1	93.6