

S2 Table . Primers used for QRT-PCR to detect the mRNA of peroxidases of *Oryza meyeriana*

Spot No	Protein name	RGAP Accession number	Sequence (5'-3')	Product size (bp)
M52	Peroxidase 4	Os07g0677100-01	ACGCCCTACTCCTTCGACAACGC	168
			AAGCCTGATCTGTCCCTGGCTCCC	
M111	Peroxidase 5	Os07g0665200-01	GACCTCTGTGACGGGAAGTGTCT	189
			CATGGCGGTTCTCATCTTGTGG	
Mock	Actin	Os11g06390.1	GAGTATGATGAGTCGGGTCCAG	125
			ACACCAACAATCCCAAACAGAG'	