

Eye exams on your mobile



Visual Acuity

Test distance vision



SightSim™

Simulate visual acuity



Tutorial

Learn how to use Peek

Units


- logMAR (1.0)
- Snellen - metric (6/60)
- Snellen - imperial (20/200)

Instructions


- Show full instructions on next test

Feedback


- Send anonymous feedback




Ask the client to show you which way the legs or the E are pointing.



When the test is finished, the phone will vibrate and you will hear an alert. The results will be shown on the screen.



Hold the phone at the clients eye height
at a distance of 2 meters.



The tested vision level is 6/60. We can show an approximation of how this is affecting your view of the world with SightSim.

