

Table S1. Antibody Source and Dilution used for Immunohistochemistry

Marker	Supplier	Clone	Dilution
Cytokeratins	Dako	AE1/AE3	1:100
CD10	Leica microsystems	56C6	1:200
VEGF	Santa Cruz Biotechnology	A20	1:100
VEGF-R2	Santa Cruz Biotechnology	A-3	1:50
B-catenin	BD	14/Beta-Catenin	1:100
E-cadherin	BD	36/E-Cadherin	1:120
Vimentin	Dako	V9	1:150
c-kit	Dako	Polyclonal	1:300

Secondary antibody: horseradish peroxidase (HRP)-conjugated polymer anti- mouse and anti-rabbit Universal DAB Detection Kit and ultraViews Universal Alkaline Phosphatase Red Detection Kit, Ventana Medical Systems Inc., Tucson, AZ, U.S.A.