

SUPPLEMENTARY TABLE 1. Flow cytometry antibodies utilized in the AIM assay for tonsils.

AIM Panel for Tonsils

Antibody	Clone	Vendor
CD4 BV650	OKT4	Biolegend
CD19 e780	H1B19	eBioscience
CD14 e780	61D3	eBioscience
CD16 e780	eBioCB16	eBioscience
CD8a e780	RPA-T8	eBioscience
CD45RA BV570	HI100	Biolegend
CCR7 APC	G043H7	Biolegend
CD25 PeCy7	BC96	eBioscience
OX40 FITC	BER-ACT35	Biolegend
PD-L1 PE	29E.2A3	Biolegend
CXCR5 BV421	J252D4	Biolegend
PD-1 BV785	EH12.2H7	Biolegend
Live/Dead e780	-	eBioscience

AIM Panel for PBMCs

Antibody	Clone	Vendor
CD4 APCeF780	RPA-T4	eBioscience
CD3 Alexa Fluor 700	UCHT-1	BD
CD8 V500	RPA-T8	BD
CD19 V500	H1B19	BD
CD14 V500	M5E2	BD
CD25 FITC	M-A251	BD
CD45RA eF450	HI100	eBioscience
CCR7 PerCPCy5.5	6043H7	Biolegend
OX40 PECy7	BER-ACT35	Biolegend
PD-L1 PE	29E.2A3	Biolegend
CXCR5 BV605	J252D4	Biolegend
Live/Dead Aqua	-	eBioscience