

S5 Table. Primer Sequence for large sequence polymorphism LSP Assays

Region of Difference	Sequence (5'-3')
RD4-Flanking Forward	CTC GTC GAA GGC CAC TAA AG
RD4-Flanking Reverse	AAG GCG AAC AGA TTC AGC AT
RD9-Flanking Forward	ACT CCC AGC GCT CGG CGG TGA CGG TAT CGT
RD9-Flanking Reverse	ATT CCG TGG GCG CTG CGG CCA ATG TTT GTT
RD12-Flanking Forward	GCC ATC AAC GTC AAG AAC CT
RD12-Flanking Reverse	CGG CCA GGT AAC AAG GAG T
RD702-Flanking Forward	TTC CGA GGA CCC GTT GTT GAG TGC
RD702-Flanking Reverse	GGG CGG GTT GGG TTG CTG GTC
RD711-Flanking Forward	GGC CGC CCT GCT CAA GAA CCT
RD711-Flanking Reverse	CCT AGG CCG GCG ACG AAG TG

The table contains a list of primer sequences used for the LSP assay