

<b>Antibody</b>	<b>Pretreatment</b>	<b>Dilution</b>	<b>Clone</b>	<b>Supplier</b>
CK7	pH 6	1:50	OV-TL 12/30 (1)	Dako, Hamburg, Germany
CK5/6	pH 9	1:50	D5/16B4	Dako, Hamburg, Germany
p40	pH 6	1:100	polyclonal	Diagnostic BioSystems, Pleasanton, USA
CK15	pH 6	1:200	polyclonal	Sigma-Aldrich, St. Louis, USA
HSP27	pH 6	1:1000	G3.1	Abcam, Cambridge, UK