

4
5
6

Supplementary Table 2. List of antibodies/reagents used in the flow cytometry analysis

Target protein	Fluorochrome	Clone	Manufacturer	Dilution	
CD27	FITC	M-T271	Biolegend	Manufacturer recommended volume	
CD19	PE	HIB19	Biolegend		
CD38	PerCP	HIT2	Biolegend		
IgD	PE-Cy7	IA6-2	Biolegend		
CD71	APC or FITC	CY1G4	Biolegend		
CD3	Alexa Fluor 700	HIT3a	Biolegend		
CD20	APC-Cy7	2H7	Biolegend		
TLR10	PE	3C10C5	Biolegend		
CD14	Alexa Fluor 700	HCD14	Biolegend		
CD16	Alexa Fluor 700	3G8	Biolegend		
CD52	PE	HI186	Biolegend		
Pax5	Alexa Fluor 488	1H9	Biolegend		1:400
IRF4	FITC	3E4	eBioscience		1:250
IRF8	PerCP-eFluor710	V3GYWCH	eBioscience	1:500	
Streptavidin	APC		Life Technologies	1:1000	
Live/Dead	Amcyan/ Aqua		Invitrogen	1:100	

7
8