

S3 Table. Antibody combinations for cell sorting

		Antibody	Sorted Cell			
			TTreg (Total T-reg)	MemTeff (Mem central)	62L-MTeff (Mem effector)	NK-CD56 (CD56 bright)
Fluorochrome	APC	CD25 (2 clones)	4 x test volume			
	Pacific Blue	CD45RA	4 x test volume			
	PE	CD56	4 x test volume			
	PE/Cy7	CD127	4 x test volume			
	BV605/eFluor605	CD62L	4 x test volume			
	FITC	TCRab	2.5 x test volume			
	AF700	CD4	4 x test volume			