

S5 Table. MIC values from CellDirector 3D, Etest and macrodilution for VSSA and hVISA.

MIC values for VSSA and hVISA strains determined by CellDirector 3D, Etest and

Macrodilution from pure cultures and from spiked blood culture bottles.

	Replicate #	Etest	Macrodilution	CellDirector 3D	
				without blood	with blood
VSSA	1	1.5	0.5	1.3	1.4
	2	1.5	1	0.86	0.93
	3	1.5	1	1.1	1.2
hVISA	1	2	1	2.9	2.5
	2	2	2	2.5	2.3
	3	2	2	1.9	1.7