

**Table S3.** List of primers used for quantitative real-time PCR reactions.

<b>Primer</b>	<b>Sequence (5'→3')</b>	<b>Amplicon size (bp)</b>
<b>18S rRNA Fwd</b>	GTA ACC CGT TGA ACC CCA TT	151
<b>18S rRNA Rev</b>	CCA TCC AAT CGG TAG TAG CG	
<b>CXCR4 Fwd</b>	TCT ACA GCA GTG TCC TCA TCC TGG CC	150
<b>CXCR4 Rev</b>	GAA TAG TCA GCA GGA GGG CAG GGA T	
<b>FOXA2 Fwd</b>	CGA GTT GAG CCT GTG AGG CCT GG	150
<b>FOXA2 Rev</b>	ACT CGG GCA ACA TGT TCG AGA ACG G	
<b>SOX17 Fwd</b>	ACC TAC AGC TAC GCG CAG GTC TCG	151
<b>SOX17 Rev</b>	TCG CGC CGT AGT ACA CGT GAA GGG	