
Supplemental Table 1. Sequences of Forward and Reverse Primers and Oligo Probes Used in TaqMan RT-PCR

Gene	Oligonucleotide	Sequence
<i>GaRDL1</i>	Forward	5'-ATACGTTTTTCATCTGACAAGTTGC-3'
	Reverse	5'-CAGCTGCTATTGTATACTTTTGCA-3'
	Probe	5'-AGTGCGAACAGCCAGCGATTGAAG-3'
<i>GaMYB2</i>	Forward	5'-CAGCATTGATGTGGACGAATTCT-3'
	Reverse	5'-CCAAGGAAACATGTTAATTACAAACT-3'
	Probe	5'-ACGGAAGGTTGCTTTGGTTTGGATT-3'
<i>GaHOX3</i>	Forward	5'-ATTTCTTCAAAGATGAAAGAACTAGA-3'
	Reverse	5'-GGCTCGTAAAACGGATATACAGT-3'
	Probe	5'-CAACGGCAACGCAGTCCAAGAAGT-3'
<i>Histone 3</i>	Forward	5'-CGGTGGTGTGAAGAAGCCTCAT-3'
	Reverse	5'-AATTTACGAACAAGCCTCTGGAA-3'
	Probe	5'-CGAAGAGCAACAGTTCCAGGACGGTAT-3'
