

1 **Supplemental Materials and Methods**

2 **Antibodies:**

3 For western blotting, anti-P54nrb (sc-376865), anti-hnRNPk (sc-28380), and anti-GAPDH (sc-
4 32233) antibodies were purchased from Santa Cruz Biotechnology. Anti-FUS (ab84078), anti-
5 PSPC1 (ab104238), anti-CIRBP (ab191885), anti-RUNX3 (ab135248), and anti-p53 (ab28) were
6 purchased from Abcam. Anti-c-Myc (9402) antibody was purchased from cell signaling. Anti-
7 PARP (4338-MC-50) antibody was purchased from TREVIGEN.

8

9 For IF, anti-Fibrillarin antibody (ab5821), anti-RPL5 (ab186857), and anti-BrdU (ab6326)
10 antibodies were purchased from Abcam. Anti-P54nrb antibody (sc-376865) was purchased from
11 Santa Cruz Biotechnology. Anti-PSF antibody (p2860) was purchased from Sigma-Aldrich.

12

13 For RNA-IP, anti-P54nrb antibody (05-950) was purchased from Millipore, and anti-PSF
14 antibody (p2860) was purchased from Sigma-Aldrich.

15

16 **Primer and probe sequences:**

17 The primers for NEAT1 RNA-FISH probe sub-cloning are: 5'-hNEAT1-Forward: 5' - TAG TTG
18 TGG GGG AGG AAG TG -3' and 5'-hNEAT1-Reverse: 5' - TAA TAC GAC TCA CTA TAG
19 GGT GGC ATG GAC AAG TTG AAG A -3'. Underlined italics indicate T7 promoter sequence.

20

21 The primer-probe sets for qRT-PCR are: NEAT1-Forward: 5' - TGG TAA GCC CGG GAC AGT
22 AA -3'; NEAT1-Reverse: 5' - CAG CGG GAA GGC CTC TCT -3'; NEAT1-Probe: 5' - CCG
23 AGT GGC TGT TGG AGT CGG TAT TG -3'; c-Myc-Forward: 5' - TGG TCT TCC CCT ACC

24 CTC TCA -3'; c-Myc-Reverse: 5' - AGA ATC CGA GGA CGG AGA GAA -3'; c-Myc-Probe:
25 5' - CGA CAG CAG CTC GCC CAA GTC C -3'; Pre-NCL1-Forward: 5' - CTC TGT CAC TGG
26 TAT CTT TTC CC -3'; pre-NCL1-Reverse: 5' - CAA AAC CAA ACC TAG AAC ACC AAA
27 TG -3'; Pre-NCL1-Probe: 5' - CAA GGC TAC TTT CTG TGG GAT GGC T -3'; LDLr-Forward:
28 5' - CAT TCA CCA AAT GAT GCC ACT T -3'; LDLr-Reverse: 5' - CGG GTG TCT CAG GCA
29 CTT AAT AA -3'; LDLr-Probe: 5' - AGA GCC TGA GTC ACC GGT CAC CCT TAA -3'; 28S-
30 Forward: 5' - CAG GTC TCC AAG GTG AAC AG -3', 28s-Reverse: 5' - CTT AGA GCC AAT
31 CCT TAT CCC G -3', 28S-Probe: 5' - TCC CTT ACC TAC ATT GTT CCA ACA TGC C -3';
32 47S rRNA precursor-Forward: 5' - TGT CAG GCG TTC TCG TCT C -3'; and 47S rRNA
33 precursor-Reverse: 5' -AGCACGACGTCACCACATC -3'. Primer probe sets for p21^{Waf1/Cip1} and
34 Gadd45a were purchased from Thermo Fisher Scientific.

35

36 **siRNAs:**

37 PSF siRNA-sense: 5'-CAG UCA UUG UGG AAC CAC UUG AAC A-5', and anti-sense 5'-
38 UGU UCA AGU GGU UCC ACA AUG ACU G-3'; FUS siRNA-sense: 5' - CGG GAC AGC
39 CCA UGA UUA AUU UGU A -5', and anti-sense 5'-UAC AAA UUA AUC AUG GGC UGU
40 CCC G -3'; RUNX3 siRNA-sense: 5' - CCC UGA CCA UCA CUG UGU UCA CCA A -5', and
41 anti-sense 5' - UUG GUG AAC ACA GUG AUG GUC AGG G -3'; CIRBP siRNA-sense: 5'-
42 GGA GGC UCC AGA GAC UAC UAU AGC A-5', and anti-sense 5'-UGC UAU AGU AGU
43 CUC UGG AGC CUC C -3'. siRNA for P54nrb (s9612 and s9614) and PSPC1 (s30594) were
44 purchased from Life Technologies.

45

46 **ASOs:**

47 ASOs used in this study all are 5-10-5 gapmer design with 10 deoxyribonucleotides in the middle
48 flanked at both ends by five 2'-O-methoxyethyl (2'-MOE) modified ribonucleotides, and is fully
49 phosphorothioate (PS) modified. NEAT1 ASO-1: 5'-ATG GGC TCT GGA ACA AGC AT-3';
50 NEAT1 ASO-2: 5'-CCA TAA GCT CAT CTG CAA AC-3'; Control ASO-1 (non-target): 5'-
51 CCT TCC CTG AAG GTT CCT CC-3'; Control ASO-2 (Pten): 5'-CTG CTA GCC TCT GGA
52 TTT GA-3'; Control ASO-3 (Malat1): 5'-TGC CTT TAG GAT TCT AGA CA-3'.

53

54 **RNA-FISH and IF:**

55 Cells seeded on glass-bottom culture dishes (MatTek) were fixed at room temperature for 30 min
56 with 4% formaldehyde in PBS and permeabilized for 5 min with 0.1% Triton X-100 in PBS. For
57 NEAT1 FISH, fluorescently labeled RNA probes were generated using FISH Tag™ RNA
58 Multicolor Kit (Thermo Fisher Scientific) according to the manufacturer's instructions. After
59 washing twice with 1×PBS and once with 2×saline-sodium citrate (SSC), permeabilized cells
60 were incubated for overnight at 55 °C with RNA probes in hybridization solution (50%
61 formamide, 2×SSC, 5% dextran sulfate, 0.5 M EDTA, 0.01% Tween-20, 100 µg/mL yeast tRNA
62 , and 1X Denhardt's salt). Cells were washed twice in wash buffer (0.2× SSC and 50%
63 formamide) at 55 °C for 30 min each, and mounted using Prolong Gold anti-fade reagent with
64 DAPI (Molecular Probes). For indirect IF staining, fixed and permeabilized cells were blocked in
65 blocking buffer (1 mg/ml BSA in PBS) at room temperature for 30 mins, and incubated with
66 primary antibodies in blocking buffer at room temperature for 2 hrs. After washing three times in
67 washing buffer (0.1% Tween-20 in PBS) for 5 mins each, cells were incubated with secondary
68 antibodies in blocking buffer at room temperature for 1 hrs. Finally, cells were washed three
69 times in washing buffer for 5 mins each and mounted using Prolong Gold anti-fade reagent with

70 DAPI (Molecular Probes). Images were generated by Confocal Laser Scanning Biological
71 Microscopy FV1000 Fluoview (Olympus) and analyzed using Fluoview Ver. 2.0b Viewer
72 (Olympus).

73

74 **WST-8 assay**

75 WST-8 assay was performed using Cell Counting Kit – 8 (Sigma) according to manufacturer's
76 instructions.

77