



Supplementary Figure 2: Repetition Experiments of [35S]-PfACTI added to *P. falciparum* cell extract. Repeated native PAGE analysis of [35S]-PfACTI translated in (A) rabbit reticulocyte (two additional replicates) or (B) *E. coli* lysate (three additional replicates) for 40 minutes (Pre), added to extract preparation buffer (mock) or *P. falciparum* cell extract (Pf lysate) at a 1:10 ratio of [35S]-PfACTI:*P. falciparum* cell extract for 15 min. Square bracket indicates possible complex formation in *P. falciparum* lysate. Arrow indicates monomeric [35S]-PfACTI.