





Supplementary Data



SUPPLEMENTARY FIG. S1. GC-MS chromatograms and mass spectra of *cis*-ocimene and *trans*-ocimene in reference mix and representative Jack Herer sample. **(A)** Ocimene peaks in CannTerpMix 1 (*cis*-ocimene 5.78 min, *trans*-ocimene 6.06 min). **(B)** Representative Jack Herer sample with α -phellandrene 4.94 min, 3-carene 5.06 min, α -terpinene, 5.24 min, and *trans*-ocimene 6.06 min. **(C)** Mass spectrum of *cis*-ocimene standard. **(D)** Mass spectrum of *trans*-ocimene standard. **(E)** Mass spectrum of *trans*-ocimene in Jack Herer sample. GC-MS, gas chromatography–mass spectrometry.





SUPPLEMENTARY FIG. S1. (Continued).

