

Supplemental Table 6. Primers used for gene analysis by quantitative Real-Time PCR.

Gene	Sequence
AREG	1: 5'-CCTCCTTCTTTCTTCTGTTTCTCC-3' 2: 5'-GTCACTATCTTTGTCTCTGCCA-3'
ERBB2	1: 5'-ACAGACCAGAGTGCAGGAT-3' 2: 5'-CCTCATCACCTACAACACAGAC-3'
THBS1	1: 5'-ACCACGTTGCTGAATTCAT-3' 2: 5'-AGTTCCTGATGGTGAATGCTG-3'
GAPDH	1: 5'-TGTAGTTGAGGTCAATGAAGGG-3' 2: 5'-ACATCGCTCAGACACCATG-3'