## Supplementary information, Figure S1




Figure S1 Size exclusion chromatogram and SDS-PAGE profiles of AcrF1/2-Csy32nt complex and AcrF1/2-Csy20nt complex. The samples were analyzed by a Superdex 200 10/600 pg size-exclusion column and eluted fractions were analyzed by SDS-PAGE. The elution volume of each sample is labeled aside the peak and all protein component bands are indicated accordingly in the SDS-PAGE photos. (A) and (B) represent the complex with 32nt and 20nt spacers, respectively.