

Table S3. RT-qPCR primers.

This Table relates to Figures 2, 4, 6, and the Methods.

PRIMERS	SOURCE
q-PCR primers for Terc, F:CCGCTGTTTTTCTCGCTGAC; R: GTTTTTGAGGCTCGGGAACG	Integrated DNA Technologies
q-PCR primers for GAPDH, F:CGTCCCGTAGACAAAATGGT; R:TTGATGGCAACAATCTCCAC	Integrated DNA Technologies
q-PCR primers for Erdr1, F:CACAGTGATGTCACCCACGA; R:GTGAGAATCGCTCCGTCCTG	Integrated DNA Technologies
q-PCR primers for Mid1 exon-8-9, F:GATCGTGCCCAACATCAAGC; R:TCGACAGTCAGGTTGTCGTG	Integrated DNA Technologies
q-PCR primers for Asmt, F: GAAGTGGGACAGGAAGTGAG; R: CGGGAACAGGAAGTGGC	Integrated DNA Technologies
q-PCR primers for Wls, F: CCAGTCTAATGGTGACCTGGG; R: TGAGAGTCAGCATGCACCAG	Integrated DNA Technologies
q-PCR primers for Tmx3, F: TACCGAGGACCACGGACTAA; R: AATACACGGTGCCTCTTCCG	Integrated DNA Technologies

q-PCR primers for Lphn2, F: ACCTCCCTGATGCCTTCAAA; R: GCTCCATGTAAGGGACACAT	Integrated DNA Technologies
q-PCR primers for Nfib, F: GGGTGGAAACCGGTGAATCTT; R: CGGAAAGGATCGGAAGGACC	Integrated DNA Technologies
q-PCR primers for Fyco1, F: ATGTCCAAAATGCCTCCGGT; R: AAGGCTCAGGAGAACAGCAC	Integrated DNA Technologies