


Phospho DNAPKLS - second gel

AJS08 12.9.2015

1°: 1:1,000 α -DNAPKLS
2°: 1:5,000 α -mouse




1. Φ TX
2. MDV3100
3. 4 Gy 1 min
4. 4 Gy 1 min + MDV
5. 4 Gy 2 min
6. 4 Gy 2 min + MDV
7. 4 Gy 5 min
8. 4 Gy 5 min + MDV
9. 4 Gy 10 min
10. 4 Gy 10 min + MDV
11. 4 Gy 30 min
12. 4 Gy 30 min + MDV

1°: 1:5,000 α -GAPDH
2°: 1:5,000 α -Rabbit
30 second exp.

Total DNAPKLS AJS08

1°: 1:1,000 α -DNAPKLS

2°: 1:5,000 α -mouse
12-9-15



1. Φ TX
2. MDV
3. 4 Gy 1
4. 4 Gy 1 + MDV
5. 4 Gy 2
6. 4 Gy 2 + MDV
7. 4 Gy 5 min
8. 4 Gy 5 min + MDV
9. 4 Gy 10 min
10. 4 Gy 10 min + MDV
11. 4 Gy 30 min
12. 4 Gy 30 min + MDV