

1 **Supplementary information**

2

3 **Biochips for Direct Detection and Identification of Bacteria in Blood Culture-**

4 **Like Conditions**

5

6

7

8 V. Templier¹, T. Livache¹, S. Boisset^{2,3}, M. Maurin^{2,3}, S. Slimani¹, R. Mathey¹, Y. Roupiez^{1*}

9 ¹ Univ. Grenoble Alpes, CNRS, CEA, INAC, SYMMES, F-38000 Grenoble

10 ² Laboratoire de Bactériologie, Centre Hospitalier Universitaire Grenoble Alpes, Univ. Grenoble Alpes, F-38000
11 Grenoble, France.

12 ³ Equipe Therex, TIMC-IMAG, CNRS UMR 5525, Univ. Grenoble Alpes, F-38000 Grenoble, France.


13

14 Corresponding author: Y. Roupiez, yoann.roupiez@cea.fr


15

16

17




20 **Figure S1:** SPRi data analysis for an initial bacterial load of $124 \pm 16 \text{ CFU.mL}^{-1}$ of *Salmonella*
 21 Enteritidis in 1 mL culture media containing human blood (10:1 ratio). **A)** Raw SPRi data. **B)**
 22 Smoothed SPRi data. **C)** First derivatives calculated from SPRi data. **D)** Second derivatives
 23 calculated from SPRi data. Yellow curves represent the averaged signal of the control plus or minus
 24 five standard deviations.



25

26 **Figure S2:** SPRi data treatment for the sterile control in experimental conditions corresponding to
 27 the Figure 2-D (culture medium and human blood in a 10:1 ratio). **A)** Raw SPRi data. **B)** Smoothed
 28 SPRi data. **C)** First SPRi derivative. **D)** Second SPRi derivative. Yellow curves represent the
 29 averaged SPR signal of the control plus or minus five standard deviations.



30

31 **Figure S3:** SPRi data treatment for the detection of *Salmonella* Enteritidis in blood culture-like
 32 conditions (*i.e.* 6 mL of blood diluted in 24 mL of culture medium, giving a 1:5 blood-to-broth ratio
 33 and a final volume of 30 mL). The initial bacterial load is 43 ± 3 CFU.mL⁻¹ representing
 34 approximately 260 bacteria per blood milliliter. **A)** Smoothed SPRi signals. **B)** Averaged SPRi data
 35 for the triplicates. **C)** First SPRi derivatives. **D)** Second SPRi derivatives. Yellow curves represent
 36 the averaged SPRi signal of the control plus or minus five standard deviations.