

1 **Supplementary materials**

2

3 **Accurate detection of methicillin-resistant *Staphylococcus aureus* in mixtures**
4 **utilizing single bacterial duplex droplet digital PCR**

5 Jun Luo^{1,2}, Junhua Li¹, Hang Yang¹, Junping Yu¹, Hongping Wei^{1#}

6 1.Key Laboratory of Special Pathogens and Biosafety, Wuhan Institute of Virology,
7 Chinese Academy of Sciences, Wuhan 430071, China

8 2. University of Chinese Academy of Sciences, Beijing 100039, China

9 #To whom correspondence should be addressed.

10 Email: hpwei@wh.iov.cn; Tel/Fax: +86 27 51861076

11

12

13

14

15

16

17

18

19

20

21

22

23 **Table S1** Bacteria strains used in this study.

Organism group	Strain	SCC <i>mec</i> type
MSSA	91118	—
MRSA	N315	II
	ZX3	II
	ZX108	I
	FY16	III
	WH70	V
	YN22	IV
	KQ6	V
	FY17	III
	ZX54	IV
MR-CoNS (<i>S. haemolyticus</i>)	WH01	
Control bacteria	<i>E. coil</i>	—

24 —: not applicable (no *mecA* gene)

25

26

27

28

29

30

31

32 **Table S2** The MRSA Index Ratio (MIR) of MRSA N315 under different concentrations.

33

Sample N315	Average <i>mecA</i> - positive droplets No	Average <i>Nuc</i> - positive droplets No	Average duplex positive droplets No	MRSA Index Ratio
1x10 ⁶ CFU/mL	7988	7562	6806	90±1.175
1x10 ⁵ CFU/mL	629	607	548	90.26±0.13
1x10 ⁴ CFU/mL	93	84	77	91.67±0.705
1x10 ³ CFU/mL	4	4	4	100

34

35

36

37

38

39

40

41

42

43

44

45

46

47

48

49 **Table S3** The MRSA Index Ratio (MIR) of mixtures of MSSA and MR-CoNS under different

50 concentrations.

51

Mixture No	Average positive droplets	<i>mecA</i> - No	Average <i>Nuc</i> - droplets No	positive	Average positive droplets No	duplex droplets	MRSA Ratio	Index
1	1 x 10 ⁶ CFU/mL WH01 + 1 x 10 ⁶ CFU/mL 91118							
	2919		5141		952		29.48±4.07	
2	1 x 10 ⁵ CFU/mL WH01+ 1 x 10 ⁵ CFU/mL 91118							
	240		679		16		6.75±0.95	
3	1 x 10 ⁴ CFU/mL WH01+ 1 x 10 ⁴ CFU/mL 91118							
	45		77		0		0	
4	1 x 10 ³ CFU/mL WH01+ 1 x 10 ³ CFU/mL 91118							
	1		2		0		0	
5	1 x 10 ⁶ CFU/mL WH01+ 1 x 10 ⁵ CFU/mL 91118							
	1696		253		52		20.74±1.56	
6	1 x 10 ⁶ CFU/mL WH01+ 1 x 10 ⁴ CFU/mL 91118							
	2575		76		14		18.24±1.24	
7	1 x 10 ⁶ CFU/mL WH01+ 1 x 10 ³ CFU/mL 91118							
	2511		6		0		0	
8	1 x 10 ⁵ CFU/mL WH01+ 1 x 10 ⁴ CFU/mL 91118							
	322		144		3		2.1±0.45	
9	1 x 10 ⁵ CFU/mL WH01+ 1 x 10 ³ CFU/mL 91118							
	349		7		0		0	
10	1 x 10 ⁴ CFU/mL WH01+ 1 x 10 ³ CFU/mL 91118							
	56		2		0		0	
11	1 x 10 ⁵ CFU/mL WH01+ 1 x 10 ⁶ CFU/mL 91118							
	661		4140		142		21.48±1.24	
12	1 x 10 ⁴ CFU/mL WH01+ 1 x 10 ⁶ CFU/mL 91118							
	85		4160		19		22.29±1.13	
13	1 x 10 ³ CFU/mL WH01+ 1 x 10 ⁶ CFU/mL 91118							
	3		1920		0		0	
14	1 x 10 ⁴ CFU/mL WH01+ 1 x 10 ⁵ CFU/mL 91118							
	68		557		2		2.9±0.6	
15	1 x 10 ³ CFU/mL WH01+ 1 x 10 ⁵ CFU/mL 91118							
	8		1426		0		0	
16	1 x 10 ³ CFU/mL WH01+ 1 x 10 ⁴ CFU/mL 91118							
	4		20		0		0	

52

53

54 **Table S4** The detection results of the duplex ddPCR and the duplex qPCR for 104 nasal specimens.

Sample	qPCR (Ct)		Result	ddPCR (droplet)			MRSA ratio	Result	Culture result
	<i>mecA</i>	<i>nuc</i>		<i>mecA</i>	<i>nuc</i>	Dual positive			
1	-	-	-	9	0	0	0	-	-
2	36.71	-	-	34	1	0	0	-	-
3	32.91	-	-	1410	8	6	75	+	+
4	34.02	-	-	30	4	0	0	-	-
5	31.21	-	-	209	2	0	0	-	-
6	31.29	-	-	127	2	0	0	-	-
7	33.41	33.59	+	160	84	82	97.6	+	+
8	37.16	-	-	14	2	0	0	-	-
9	-	-	-	3	2	0	0	-	-
10	-	-	-	6	4	0	0	-	-
11	-	-	-	2	2	2	100	+	+
12	-	-	-	0	0	0	0	-	-
13	34.72	35.91	+	13	5	0	0	-	-
14	29.98	32.37	+	3052	1542	1483	96.2	+	+
15	32.01	32.51	+	463	450	379	84.22	+	+
16	37.39	-	-	0	0	0	0	-	-
17	-	-	-	2	0	0	0	-	-
18	-	-	-	3	3	0	0	-	-
19	-	-	-	1	1	1	100	+	+
20	37.67	-	-	0	0	0	0	-	-
21	-	-	-	0	0	0	0	-	-
22	-	-	-	1	0	0	0	-	-
23	-	-	-	1	0	0	0	-	-
24	35.17	-	-	195	1	1	100	+	+
25	35.6	37.83	+	10	8	0	0	-	-
26	-	-	-	0	3	0	0	-	-
27	33.77	-	-	87	0	0	-	-	-
28	-	-	-	1	1	0	-	-	-
29	34.13	-	-	69	1	0	-	-	-
30	36.01	-	-	9	0	0	-	-	-
31	-	-	-	0	2	0	-	-	-
32	37.12	-	-	7	3	0	-	-	-
33	34.35	-	-	117	3	2	66.67	+	+
34	-	-	-	21	1	0	-	-	-
35	-	-	-	0	0	0	-	-	-
36	35.74	-	-	14	5	0	-	-	-
37	30.32	-	-	292	1	0	-	-	-
38	-	-	-	0	6	0	-	-	-

39	-	-	-	3	0	0	-	-	-
40	34.01	-	-	15	0	0	-	-	-
41	30.95	37.74	+	51	9	8	88.89	+	+
42	33.94	-	-	21	3	0	-	-	-
43	32.18	-	-	33	0	0	-	-	-
44	-	-	-	3	0	0	-	-	-
45	31.14	-	-	55	2	2	100	+	+
46	-	-	-	0	0	0	-	-	-
47	-	-	-	0	0	0	-	-	-
48	31.62	-	-	26	1	0	-	-	-
49	30.38	-	-	107	3	3	100	+	+
50	36.89	-	-	4	0	0	-	-	-
51	35.14	-	-	8	2	0	-	-	-
52	33.21	-	-	7	0	0	-	-	-
53	-	-	-	2	1	0	-	-	-
54	36.56	-	-	6	2	0	-	-	-
55	-	-	-	3	0	0	-	-	-
56	-	-	-	0	0	0	-	-	-
57	-	-	-	1	0	0	-	-	-
58	-	-	-	4	3	0	-	-	-
59	-	-	-	1	3	0	-	-	-
60	-	-	-	0	0	0	-	-	-
61	-	-	-	10	0	0	-	-	-
62	-	-	-	4	0	0	-	-	-
63	32.96	-	-	33	1	0	-	-	-
64	-	-	-	1	1	0	-	-	-
65	-	-	-	0	1	0	-	-	-
66	31.71	-	-	238	1	0	-	-	-
67	-	-	-	1	1	0	-	-	-
68*	24.6	-	-	10374	10	10	.*	-	-
69	37.95	-	-	3	0	0	-	-	-
70	32.12	-	-	143	0	0	-	-	-
71	-	-	-	8	0	0	-	-	-
72	36.14	-	-	18	4	1	25	+	+
73	-	-	-	0	4	0	-	-	-
74	31.48	35.41	+	231	60	17	28.33	+	+
75	-	-	-	1	4	0	-	-	-
76	34.03	37.99	+	31	7	1	14.28	+	+
77	-	-	-	1	2	0	-	-	-
78	29.56	-	-	158	2	0	-	-	-
79	31.12	-	-	59	0	0	-	-	-
80	-	-	-	0	0	0	-	-	-
81	-	-	-	0	0	0	-	-	-
82	-	-	-	4	3	3	100	+	+

83	36.21	36.65	+	11	10	7	70	+	+
84	-	-	-	2	1	0	-	-	-
85	35.1	-	-	29	0	0	-	-	-
86	-	-	-	6	0	0	-	-	-
87	35.62	-	-	11	0	0	-	-	-
88	-	-	-	3	0	0	-	-	-
89	-	-	-	1	4	0	-	-	-
90	37.19	-	-	8	0	0	-	-	-
91	31.91	-	-	157	0	0	-	-	-
92	-	-	-	5	0	0	-	-	-
93	-	-	-	5	7	3	60	+	+
94	31.56	-	-	230	0	0	-	-	-
95	32.79	-	-	151	0	0	-	-	-
96	32.61	-	-	145	1	1	100	+	+
97	-	-	-	1	1	0	-	-	-
98	-	-	-	1	0	0	-	-	-
99	-	-	-	1	1	0	-	-	-
100	37.23	-	-	11	1	0	-	-	-
101	-	-	-	1	1	0	-	-	-
102	31.98	-	-	125	5	0	-	-	-
103	-	-	-	2	0	0	-	-	-
104	33.1	-	-	101	1	0	-	-	-

55 **Note:** “-“ means no signal (for qPCR) or MRSA negative result. “+“ means MRSA positive result.

56 ***: Specimen 68 was diluted 10 times for determining the MIR in the ddPCR assay.**

57

58

59

60

61

62

63

64

65

66

67

68

69

70


71

72


73

74 Fig.S1

A



B



75

76 Fig.S1 One-dimensional (1D) plots of the duplex ddPCR assay under different annealing

77 temperatures (54–60°C). MRSA strain N315 (1×10^7 CFU/mL) is used as the template (*A. mecA*

78 gene and *B. nuc* gene). Each point represents a droplet with a given fluorescence level and the

79 droplet colors indicate which target was amplified. The x-axis shows various temperature, and the

80 y-axis represents the fluorescence amplitude.

81

82

83

84

85

86

87

88

89

90


91

92

93

94 Fig.S2

95



96

97 Fig. S2 The activity changes of ClyH after cooling on ice. The concentration of ClyH:

98 10 µg/mL.

99

100

101

102

103

104

105

106

107

108

109

110

111

112

113

114

115


116

117

118

119

120
121
122



123


124 **Fig. S3** Two-dimensional (2D) plots of the duplex ddPCR for detection of the mixture of
125 MSSA(91118) and MR-CoNS under different concentrations (Sample No 5 - No 8 are as shown in
126 Table S3). The x-axis shows the fluorescence amplitude corresponding to the HEX fluorophore
127 (*nuc*), and the y-axis represents the fluorescence amplitude corresponding to the FAM fluorophore
128 (*mecA*). Each point represents a droplet with a given fluorescence level and the droplet colors
129 indicate which target was amplified (blue: *mecA* positive; green: *nuc* positive; black: negative,
130 Orange: both *nuc* and *mecA* positive).

131

132

133

134



135

136 **Fig. S4** Two-dimensional (2D) plots of the duplex ddPCR for detection of the mixture of
 137 MSSA(91118) and MR-CoNS under different concentrations (Sample No 9 - No 12 are as shown
 138 in Table S3). The x-axis shows the fluorescence amplitude corresponding to the HEX fluorophore
 139 (*nuc*), and the y-axis represents the fluorescence amplitude corresponding to the FAM fluorophore
 140 (*mecA*). Each point represents a droplet with a given fluorescence level and the droplet colors
 141 indicate which target was amplified (blue: *mecA* positive; green: *nuc* positive; black: negative,
 142 Orange: both *nuc* and *mecA* positive).

143

144


145

146

147

148

149



150

151 **Fig. S5** Two-dimensional (2D) plots of the duplex ddPCR for detection of the mixture of
 152 MSSA(91118) and MR-CoNS under different concentrations (Sample No 13 - No 16 are as shown
 153 in Table S3). The x-axis shows the fluorescence amplitude corresponding to the HEX fluorophore
 154 (*nuc*), and the y-axis represents the fluorescence amplitude corresponding to the FAM fluorophore
 155 (*mecA*). Each point represents a droplet with a given fluorescence level and the droplet colors
 156 indicate which target was amplified (blue: *mecA* positive; green: *nuc* positive; black: negative,
 157 Orange: both *nuc* and *mecA* positive).

158

159

160