

**S2 Table** Intra- and interday precision of the 6-plex assay (%RSD) for fortified samples (in dark ale<sup>a</sup>)

<b>Precision</b>	<b>n</b>	<b>AFB<sub>1</sub></b> (0.5 µg/L)	<b>DON</b> (100 µg/L)	<b>FB<sub>1</sub></b> (100 µg/L)	<b>OTA</b> (2 µg/L)	<b>T-2</b> (50 µg/L)	<b>ZEN</b> (50 µg/L)
Intraday 1	11	7.5	4.3	6.7	3.7	10.9	5.9
Intraday 2	6	3.6	8.1	8.3	10.6	16.7	5.8
Interday	17	6.5	6.5	7.7	7.4	18.4	8.5

<sup>a</sup> n = the total number of fortified dark ale samples