

**S4 Table** Intra- and interday determinations of mycotoxin concentrations in a multi-fortified dark ale using the 6-plex assay<sup>a</sup>

Mycotoxin	Mycotoxin concentration (µg/L)						
	Spike	Intraday 1 (n=11)		Intraday 2 (n=6)		Interday (n=17)	
		Average	Median	Average	Median	Average	Median
AFB <sub>1</sub>	0.5	0.7	0.6	0.7	0.7	0.7	0.6
DON	100	126	122	116	125	122	124
FB <sub>1</sub>	100	286	275	385	414	321	282
OTA	2	3.3	3.3	2.9	2.9	3.2	3.2
T-2	50	48	48	54	54	50	48
ZEN	50	61	58	69	69	64	64

<sup>a</sup> n = the total number of fortified dark ale samples