



SUPPLEMENTARY FIG. S3. SEC analysis of Dox-LNP. DSPE-PEG₂₀₀₀-COOH was conjugated to fluorescent Dox. The SEC chromatogram of the Dox-LNP shows the Absorbance₂₁₄ (*solid line*) coregisters with the fluorescent 550 nm peak (*dashed line*). Dox, doxorubicin; DSPE-PEG, distearoylphosphatidyl ethanolamine-polyethylene glycol; LNP, lipid nanoparticle; SEC, size-exclusion chromatography.