

Table S1: Flow Cytometry Panels

Tube #	FITC	PerCP	APC	APC-Cy5.5/AF-700	APC-Cy7/AF-750	PE	PE-TR/PE-AF610	PE-Cy5	PE-Cy5.5	PE-Cy7
1	CD23	CD20	IgM	CD38	CD27	CD10	CD24	X	CD5	CD19
2	CD23	X	IgM	CD38	CD27	IgD	CD24	CD21	X	CD19
3	CD23	CD20	HLA-DR	CD38	CD27	X	X	CD138	X	CD19
4	IgD	9G4	IgM	CD38	CD27	X	X	X	CD5	CD19
5	CD23	X	IgM	CD38	CD27	CD267	X	CD95	X	CD19
6	NKG2D	CD45	NKp44	CD56	CD3	CD8	CD4	CD62L	X	CD16

Tubes 1-5 Live Dead Dump #1 = CD3, CD8, CD14

Tube 6 Live Dead Dump #2 = CD19, CD 14

X = blank