

Identification of novel serum autoantibodies against EID3 in non-functional pancreatic neuroendocrine tumors

SUPPLEMENTARY MATERIALS

Supplementary Table 1: List of antibodies used for immunohistochemical staining and Western blotting

Immunohistochemical staining			
Molecule	Clonality	Dilution	Manufacture
EID3	Rabbit polyclonal	1:1000	Abcam
Western blotting			
Molecule	Clonality	Dilution	Manufacture
PTEN	Mouse monoclonal	1:2000	Origene
EID3	Rabbit polyclonal	1:1000	Abcam
EHD1	Rabbit monoclonal	1:1000	Abcam
LGALS9	Rabbit polyclonal	1:1000	GeneTex
BRAP	Rabbit polyclonal	1:250	Abnova