

Set	Forward primer	Reverse primer	Purpose	Product size (bp)
1	CGTGCCTACCGAGGAGACC	GCGTAGCGCCTACATCAGCA	Circularization	460
2	CGTACCTCTGTTGGCAACCC	GTAGCTCCCGGTGGAACCAA	Circularization	373
Pos1	GCGCGAAACATAGACTGCTG	GAGAACAACCATGGACACGG	Positive control	255
Pos2	GCTTAGCCGAGTTCATTCG	CTACGGTACGTGCGCCTAGC	Positive control	292

