

Supplementary Information

One-step noninvasive prenatal testing (NIPT) for autosomal recessive homozygous point mutations using digital PCR

Mun Young Chang,^{1*} Soyeon Ahn,^{2*} Min Young Kim³, Jin Hee Han³, Hye-Rim Park³, Han Kyu Seo⁴, Jinsun Yoon⁵, Seungmin Lee³, Doo-Yi Oh³, Changsoo Kang⁶ & Byung Yoon Choi^{3,7}

¹Department of Otorhinolaryngology-Head and Neck Surgery, Chung-Ang University College of Medicine, 102 Heukseok-ro, Dongjak-gu, Seoul, 06973, Republic of Korea.

²Medical Research Collaborating Center, Seoul National University Bundang Hospital, 82 Gumi-ro 173 beon-gil, Bundang-gu, Seongnam, 13620, Republic of Korea.

³Department of Otorhinolaryngology, Seoul National University Bundang Hospital, 82 Gumi-ro 173 beon-gil, Bundang-gu, Seongnam, 13620, Republic of Korea.

⁴LAS Inc., 16 Arayuk-ro, Gimpo, 10136, Republic of Korea.

⁵Bio-Medical Science Co., Ltd., BMS Bldg., 22 Yeoksam-ro 7-gil, Gangnam-gu, Seoul 06244, Republic of Korea.

⁶Department of Biology and Research Institute of Basic Sciences, College of Natural Sciences, Sungshin Women's University, Dongseon-dong 3-ga, Seongbuk-gu, Seoul, 01133, Republic of Korea.

⁷Wide River Institute of Immunology, Seoul National University College of Medicine, 101 Dabyeonbat-gil, Hwachon-myeon, Hongcheon, 25159, Republic of Korea.

Correspondence to : Byung Yoon, Choi, MD, PhD

Department of Otorhinolaryngology, Seoul National University Bundang Hospital, 82 Gumi-ro 173 beon-gil, Bundang-gu, Seongnam 13620, Korea, Tel: +82-31-787-7406; Fax: +82-31-787-4057

Supplementary Table S1. The sequences and concentrations of primers and probes in picodroplet digital PCR.

Name	Sequence(5' to 3')	Conc. Used for RainDrop(μ M)	Anneal Temp.
c.235delC-F	CCATCTCCCACATCCGGC	0.35	60
c.235delC-R	CACGTGCATGGCCACTAG	0.35	
c.235deC-WT Probe	VIC-AGCTGCAGGGCCCAT-NFQ	0.2	
c.235delC-MT probe	FAM-ATCAGCTGCAGGGCCCAT-NFQ	0.2	

Conc, concentration; Temp, temperature.

Supplementary Table S2. The protocol of polymerase chain reaction in picodroplet digital PCR.

Step	Temp	Time	Cycles
Polymerase activation	95 °C	10 min	1
Denaturation	95 °C	15 sec	45
Annealing & Extension	*Depending on target	1 min	45
Incubation	98 °C	10 min	1
Final Hold	12 °C	Hold	

* A slow ramping speed (0.5°C) was used during cooling from the denaturation step to the annealing step.


Temp, temperature.

Supplementary Table S3. The sequences and concentrations of primers and probes in chip-based dPCR.


Name	Sequence(5' to 3')	Conc. Used for QX200 (uM)	Anneal Temp.
SLC26A4_IVS7-2A>G -F	TTGACAAACAAGGAATTATTTAAAACCAATG GA	0.35	
SLC26A4_IVS7-2A>G -R	TCCAGGTTGGCTCCATATGAAATG	0.35	62
SLC26A4_IVS7-2A>G - WT Probe	VIC-CAATTATCGTCTGAAATAA-NFQ	0.2	
SLC26A4_IVS7-2A>G - MT probe	FAM-ATCGTCCGAAATAA-NFQ	0.2	

Supplementary Table S4. The protocol of polymerase chain reaction in chip-based dPCR.


Step	Temp	Time	Cycles
Enzyme activation	95 °C	10 min	1
Denaturation	94 °C	30 sec	40
Annealing / Extension	62 °C	1 min	
Enzyme deactivation	98 °C	10 min	1
Cool down and storage	4 °C	Hold	1




Supplementary Figure S1. Pedigree of three families. (A) Sanger sequencing traces of SH197 family. Father: *GJB2* c.235delC carrier, Mother: *GJB2* c.235delC carrier, the first baby: *GJB2* c.235delC homozygote. (B) Sanger sequencing traces of SB275 family. Father: *GJB2* c.235delC carrier, Mother: *GJB2* c.235delC carrier, the first baby: *GJB2* c.235delC homozygote. (C) Sanger sequencing traces of SH162 family. Father: *SLC26A4* IVS7-2A>G carrier, Mother: *SLC26A4* IVS7-2A>G carrier, the first baby: *SLC26A4* IVS7-2A>G homozygote.




Supplementary Figure S2. Validation of probes. The target was detected after less than 40 cycles. (A) The probe for mutation in SH197 and SB275 families (*GJB2* c.235delC), (B) The probe for mutation in SH162 family (*SLC26A4* c.2168A>G).




Supplementary Figure S3. Two-dimensional histogram of the mutation (*GJB2* c.235delC) in wildtype (A,B), heterozygote (C) and homozygote (D,E) controls and maternal plasma DNA (F,G) of SB275 family.




Supplementary Figure S4. Two-dimensional histogram of the mutation (*SLC26A4* IVS7-2A>G) in wildtype (A,B), heterozygote (C) and homozygote (D,E) controls and maternal plasma DNA (F,G) of SH162 family using picodroplet digital PCR (1st trial).



Supplementary Figure S5. Two-dimensional histogram of the mutation (*SLC26A4* IVS7-2A>G) in wildtype (A), heterozygote (B,C) controls and maternal plasma DNA (D,E) of SH162 family using picodroplet digital PCR (2nd trial).



Supplementary Figure S6. Two-dimensional histogram of the mutation (*SLC26A4* IVS7-2A>G) in wildtype (A), heterozygote (B,C) controls and maternal plasma DNA (D) of SH162 family using chip-based dPCR (1st trial).



Supplementary Figure S7. Two-dimensional histogram of the mutation (*SLC26A4* IVS7-2A>G) in wildtype (A,B), heterozygote (C) and homozygote (D,E) controls and maternal plasma DNA (F,G) of SH162 family using chip-based dPCR (2nd trial).