

<b>Supplemental Table 1:</b>			
<b>Fluorophore</b>	<b>T-cell Panel</b>	<b>Clone</b>	<b>Manufacturer</b>
APC	GP66+ tetramer		
PE-Cy5	x		
AF700	CD44	IM7	eBioscience
APC-Cy7	x		
PE	Tbet	ebio 410	eBioscience
PE-Cy7	CXCR5	2g8	eBioscience
PE-CE594	b220	RA4-682	BD
PE-e610	x		
FITC	PD-1	J43	eBioscience
P5.5			
PE-e710			
BV421	CXCR3		173 BD
e450			
BV510	CD4	RM 4-5	BD
BV605	CCR7-bio	4B17	BD
BV650			
BV711	CD8	53-6.7	BD
BV786	CD62L	MEL14	BD
BUV 737	x		
BUV 395	Thy1.2	53-2.1	BD

Fluorophore	B cell Panel	Clone	Manufacturer
APC	IgM	II/41	BD
PE-Cy5	DECOY		
AF700	CD38	90	eBioscience
APC-Cy7	none		
PE	MSP tetramer		
PE-Cy7	CD73	Ty11.8	eBioscience
PE-CE594	x		
PE-e610	x		
FITC	CD80	16-10A1	Biolegend
P5.5	CD4	RM 4-5	BF
PE-e710			
BV421	x		
e450	GL7	GL7	BD
BV510	live/dead		
BV605	CD138	281-2	BD
BV650	IgD	11-26-2a	BD
BV711	CD19	1D3	BD
BV786			
BUV 737	B220	RA3-GB2	BD
BUV 395	Thy 1.2	53-21	BD
	<b>Congenic</b>		
BV650	CD45.1	A20	BD
FITC	CD45.2	104	Biolegend

<b>Fluorophore</b>	<b>Parasitemia ant</b>	<b>Clone</b>	<b>Manufacturer</b>
APC	CD45	R17217	BD
FITC	Ter119	Ter119	BD
PE	CD71	30-F11	eBioscience

<b>Fluorophore</b>	<b>Global Panel</b>	<b>Clone</b>	<b>Manufacturer</b>
APC	Siglec H	551	Biolegend
PE-Cy5	NONE		
AF700	CD4	GK1.5	eBioscience
APC-Cy7	Ly6C	HK1.4	eBioscience
PE	CD19	1D3	BD
PE-Cy7	Ly6G	1a8	biolegend
PE-CE594	CD11b	M1/70	BD
PE-e610	x		
FITC	CD317	eBio927	eBioscience
P5.5			
PE-e710			
BV421			
e450	MHCII	M5/114.15.2	eBioscience
BV510	cd8	53.617	
BV605	F4/80 bio	BM8	biolegend
BV650			
BV711	CD11c	HL3	BD
BV786	CD62L	MEL14	BD
BUV 737	B220	RA3-GB2	BD
BUV 395	Thy 1.2	53-21	BD