

Supplementary Table S1 Real-time PCR primers and probes.

Target	Accession number	Forward primer	Reverse primer	Probe
RPS10 ^a	NM_001193579.1	AATGTGCCCAACCTTCATGTC	TCCAGGCAAACCTGTTCTTCA	TGAAGCCATGCAGTCTCTCAAGTCCC
ESR1	XM_004011538.1	CTGAGCGACCTCACTTTCAC	CCAGAGTGCCTGCTCCTTT	TTTCAGGCTGTGACCAAGGGC
PGR	XM_004016530.1	GGCACTGAGTGTTGAATTTCCA	CACCATCCCTGCCAAGATC	TGATGTCTGAGGTTATCGCTGCACA
MMP26	XM_015115062.1	CCTGCAGCACTCTGGGAATC	TTCCAGGGTCGTGATACCAGTA	AGCTCCATAATGTACCCC
TIMP3	NM_001032956.1	GGACCGACATGCTCTCCAAT	CGGATGCAGGCATAGTGTGTTG	TTCGGCTACCCTGGCTACCA

^aEndogenous control primer.