

**S3 Table. Primers used for PCR analysis for RD cell line.**

Sequences for the forward (For.) and reverse (Rev.) primers for each target are shown.

Sequences are shown from 5' to 3'.

Gene	ID	Sequence	Size
<i>hDMD</i>	For.	GGACCATTTGACGTTAAGGAAAC	
1st PCR	primer		
	Rev.	AGGTCTTTGGCCAACTGCTTG	
	primer		
<i>hDMD</i>	For.	CACTCAGCCAGTGAAGAGGA	497 bp
2nd PCR	primer		(exon 51 inclusion)
	Rev.	ACTGATTCTGAATTCTTTCAATTCG	264 bp
	primer		(exon 51 skipped)
<i>hGAPDH</i>	For.	ACCACAGTCCATGCCATCAC	452 bp
	primer		
	Rev.	TCCACCACCCTGTTGCTGTA	
	primer		