

Supplementary Table S1: Elution gradient of the kynurenine pathway metabolites HPLC-MS/MS assay. 0.2% formic acid in water and acetonitrile were used as mobile phases A and B, respectively.

Step No	Time	% Mobile Phase A	% Mobile Phase B	Flow μL/min
1	0	84	15	650
2	1.2	85	15	650
3	3.2	72	28	650
4	4.5	65	35	650
5	5.2	5	95	650
6	5.5	5	95	650
7	5.6	85	15	650
8	10.5	85	15	650