

SUPPLEMENTARY TABLE S3. PCR SETTINGS FOR THE ANALYSIS OF mRNA AND miRNA LEVELS

<i>mRNA/SYBR green</i>				
<i>Hold</i>		<i>40 cycles</i>		<i>Melt curve</i>
10 min 95°C	30 s 95°C	60 s 60°C	45 s 72°C	— 60–95°C
<i>mRNA (lipid-regulated genes)/TaqMan</i>				
<i>Hold</i>	<i>Hold</i>	<i>40 cycles</i>		
2 min 50°C	10 min 95°C	15 s 95°C	60 s 60°C	
<i>miRNA/SYBR green</i>				
<i>Hold</i>	<i>Hold</i>	<i>40 cycles</i>		<i>Melt curve</i>
2 min 50°C	2 min 95°C	15 s 95°C	60 s 60°C	— 60–95°C