

Supporting information for

A Highly Sensitive Immunochromatographic Strip Test for Rapid and Quantitative Detection of Saikosaponin d

Table S1 Variations among ELISA runs for the analysis of SSd

SSd (ng mL ⁻¹)	CV%	
	ELISA ¹	ICS ²
500	3.76	2.49
1000	4.29	2.64
2000	3.47	4.38

¹For the variations, measured values indicate mean \pm SD for samples run in triplicate on three different microtiter plates. ²For the variations, six strips were tested for each concentration.