



**Figure S1. Staining controls for expression of IL-6 and IL-8 receptors.** OneComp eBeads (compensation beads from eBioscience) were used as a positive control to demonstrate successful binding of fluorochrome-conjugated antibodies in comparison to HeLa cells and OneComp eBeads incubated without fluorochrome-conjugated antibody. Left panel shows positive staining for FITC-conjugated anti-CD181 (IL-8RA chain). Middle panel shows positive staining for PE-conjugated anti-CD182 (IL-8RB chain). Right panel shows positive staining for APC-conjugated anti-CD126 (IL-6R $\alpha$  chain).