

Supplemental Data Table S4. Dilution tests for AFIAS anti-HBs test

No. Case	KFDA standard 08/026 (IU/L)	AFIAS anti-HBs (COI)*	CV %
1	25,000.0	1,521.4	3.2
2	12,500.0	1,501.3	1.9
3	6,250.0	1,307.5	1.7
4	3,125.0	1,175.2	1.9
5	1,562.5	1,039.3	7.2
6	781.3	724.4	1.9
7	390.6	389.5	6.0
8	195.3	189.5	4.2
9	97.7	102.8	7.6
10	48.8	69.3	4.0
11	24.4	28.6	8.7
12	12.2	11.3	8.9
13	6.1	5.3	8.7
14	3.1	2.5	42.3

*Mean of five measurements.

Abbreviations: AFIAS, Automated Fluorescent Immunoassay System; anti-HBs, antibody to HBsAg; KFDA, Korea Food and Drug Administration; COI, cut-off index.