

S2 Table. Sequences of primers used for real-time PCR.

Primer name	Primer Sequence	Product size (bp)	Annealing Temperature
InfA F	TTTCACCGAGGAGGGAGCA	105	56.5°C
InfA R	CCTCCGATGAGGACCCCAA		
Influenza Bs	GTCCATCAAGCTCCAGTTTT	145	
Influenza Bas	TCTTCTTACAGCTTGCTTGC		
RSVB F	GCATTAGCCAAAGCAGCAATAC	155	
RSVB R	CATCCATTAGACCATTTAATTGAGAGC		