


Supplement. Table 1. Sequences of the primers used in RT-PCR analysis.


Gene	Primer	Sequence	Accession number
COX-2	Sense primer	CAGCAAATCCTTGCTGTTCC	XM-004028064.1
	Anti-sense primer	TGGGCAAAGAATGCAAACATC	
EP1	Sense primer	ATGTCCGAAGCAAACATCAC	NM-010927.3
	Anti-sense primer	TAATGTCCAGGAAGTAGGTG	
EP2	Sense primer	ATGAGCACAGAAAGCATGATC	NM-013693.1
	Anti-sense primer	TACAGGCTTGTCACCTCGAATT	
EP3	Sense primer	CTCGTGCTGTCCGACCCATAT	NM-008361.3
	Anti-sense primer	TTGAAGACAAACCGCTTTTCCA	
EP4	Sense primer	TGGTATGTGGGCTGGCTG	NM-000958
	Anti-sense primer	GAGGACGGTGGCGAGAAT	
GAPDH	Sense primer	TTCACCACCATGGAGAAGGC	XR031141.1
	Anti-sense primer	GGCATGGACTGTGGTCATGA	

sFig. 1. Jeong *et. al.*



sFig. 2. Jeong *et. al.*

A)



B)

