


**Additional file 5.** Flow cytometric gating strategy of stimulated plasmacytoid dendritic cells and NK cells



**Additional figure S4.** Plasmacytoid dendritic cells (pDCs) and NK cells stimulated in co-cultures with RNA-containing immune complexes (RNA-IC). (A) Cells were consecutively gated for live cells, singlets and lymphocytes. (B) Cells were stained with monoclonal antibodies to blood dendritic cell antigen (BDCA) 2 and CD56 and analyzed by flow cytometry to identify pDC (red) and NK cells (blue), respectively.