

Supplementary data

S1 Table. List of primers used for RT-qPCR

Gene	Sense primer (5'-3')	Antisense primer (3'-5')
SOD-1	ACCAGTGTGCGGCCAATGATG	GCTGTACCAGTGCAGGTCCTCA
SOD-2	AAGCTGACGGCTGCATCTGTTG	CGTGCTCCCACACATCAATCCC
GPx	AGTTCGGACATCAGGAGAA	AGGGCTTCTATATCGGGTTC
CAT	AGTCTCGCCGCATCTTCAACAG	TCTGGGACTTCTGGAGCCTACG
GAPDH	GGGAGCCAAAAGGGTCATCA	TGATGGCATGGACTGTGGTC