



Figure S1. The intracellular levels of p24CA in SupT1/Ank^{GAG}1D4, and SupT1/Ank^{A32}D3 cells on day 7 post-infection. The SupT1/Ank^{GAG}1D4 and SupT1/Ank^{A32}D3 cells were collected on day 7 post-infection for this purpose. The representative 14 samples in each cluster (C1, C2, C3, and C4) and two samples (DK071 and UU80), for which Ank^{GAG}1D4 did not show the inhibition effect, were selected. Cellular proteins were extracted by RIPA buffer. The total protein lysates were normalised to 10 µg/mL to determine the intracellular levels of p24CA using Genscreen ULTRA HIV Ag-Ab Assay.

