

Supplementary Table 1. Antibodies used in a case of progressive encephalomyelitis with rigidity and myoclonus.

Antibody	Clone	Company	Code	Dilution	Detection Method
CD3	LN10	Leica Biosystems	NCL-L-CD3-565	100	ABC
CD4	4B12	Leica Biosystems	NCL-L-CD4-368	50	ABC
CD20cy	L26	Agilent Technologies	IR604	Ready-to-use	ABC
CD34	NU-4A1	Nichirei Bioscience	413111	Ready-to-use	Polymer
CD68	KP1	Agilent Technologies	IR609	2,000	Polymer
c-Kit (CD117)	rabbit poly	Agilent Technologies	A4502	200	Polymer
Desmin	D33	Monosan	MON3001	50	Polymer
DOG-1	RBT-DOG1	Bio SB	BSB6273	50	Polymer
IgG	rabbit poly	Agilent Technologies	IR512	4	Polymer
S-100A	rabbit poly	Agilent Technologies	Z0628	2,000	Polymer
α -SMA	1A4	Agilent Technologies	M0851	400	Polymer
PDGFRa	rabbit poly	Santa Cruz Biotechnology	sc-338	80	Polymer
PDGFRb	rabbit poly	Santa Cruz Biotechnology	sc-339	100	Polymer

α -SMA, Alpha smooth muscle actin; PDGFR, Platelet-derived growth factor receptor.