

Supplemental Table 3: Primers used for ChIP-qPCR

	Forward Primer	Reverse Primer
ChIP Myb (N)	Sense, 5'-CATGCACTCAGAGACAGAG-3'	antisense, 5'-GCAAATCATTCAAACGCCGC-3'
ChIP Myb (E)	Sense, 5'-GCTGGAGGCATGCTGTAAG-3'	antisense, 5'-GCCATTTTGCCTCACAACC-3'
Nfil3 (A)	Sense, 5'-ATGTCTCCATCCTTGAAGTG-3'	antisense, 5'-GCTTGCATTTTCTGTCTTTT-3'
Nfil3 (B)	Sense, 5'-TGGTTTTTGGAGACGTTTTCT-3'	antisense, 5'-GTTTAATCCCCAAAGCCTAT-3'
Nfil3 (C)	Sense, 5'-CAAGGGCTTTATCATCTGAG-3'	antisense, 5'-CCATCTGTGAGCTACTCTCC-3'
ChIP Bic (A)	Sense, 5'-GGTCCTAAAGGGAACCCACAG-3'	antisense, 5'-TGCTGTTTCCTGTTAGAGTTGG-3'
ChIP Bic (B)	Sense, 5'-AGGTCAAAGGCAACCAAAG-3'	antisense, 5'-TGGATTTCCCACCAGAATTT-3'
ChIP Wnt10a (A)	Sense, 5'-TTAATGAGAGTGGGAGCCGG-3'	antisense, 5'-CCCCAGCAGAACCCAAGG-3'
ChIP Wnt10a (B)	Sense, 5'-TGCTCAATGGTTAAGGGCAC-3'	antisense, 5'-AGAGCCAAAACATGATGCCG-3'