



**Figure S12.** HPLC/MS analyses of AFD10 crude extract. HPLC-conditions: Detection wavelength 254 nm; solvent A: H<sub>2</sub>O/0.1% Formic acid, solvent B: CH<sub>3</sub>CN/0.1% Formic acid; flow rate: 0.5 mL min<sup>-1</sup>; 0-4 min, 10% B; 4-22 min, 10-100% B; 22-27 min, 100% B; 27-29 min, 100%-10% B; 29-30 min, 10% B. (MW = Molecular Weight)