

Supplemental Table S1 – Lipid composition of SPLASH LipidoMix™, which contains a mixture of all major lipid classes found in biological samples.

Lipid Species	Concentration (µg/mL)
PC(15:0_18:1{d7})	160.7
PE (15:0_18:1{d7})	50.8
PS(15:0_18:1{d7})	203.7
PG(15:0_18:1{d7})	29.1
PI(15:0_18:1{d7})	900.8
PA(15:0_18:1{d7})	357.4
LPC(18:1{d7})	25.5
LPE(18:1{d7})	95.0
CE(18:1{d7})	356.1
MAG(18:1{d7})	199.7
DAG(15:0_18:1{d7})	9.4
TAG(15:0_18:1{d7}_15:0)	57.3
SM(18:1{d7})	30.9
Cer(15:0{d7})	19.7