

Supplementary Table 1: Successive versions of the screening tube

Version	PacB/BV421	PacO/BV510	FITC	PE	PerCPCy 5.5	PECy7	APC	APC-A750
1	CD4+IgM PacB	CD45 PacO	CD8+smlgD	CD16 +CD56	CD27	CD19+TCRgd	CD3	CD38
2	CD27 BV421	CD45 PacO	CD8+smlgD	CD16 +CD56	CD4+IgM	CD19+TCRgd	CD3	CD38
3	CD27 BV421	CD45 PacO	CD8+smlgD	CD16 +CD56	CD4+IgM	CD19+TCRgd	CD3	CD45R A
4 (Final)	CD27 BV421	CD45RA BV510	CD8+smlgD	CD16 +CD56	CD4+IgM	CD19+TCRgd	CD3	CD45