

Supplementary table 1

list of qPCR primers used in this study

<i>Designation</i>	<i>5' – 3' Forward</i>	<i>5' – 3' Reverse</i>	<i>Manufacturer</i>
SyBR Green			
TP53i9 minigene α	CTTCAGATCCGTGGGCGTGA	CTGCAGAATTCGCCCTTC TG	Eurogentec
TP53i9 minigene β	AAGAGAGCATGAAAATGGT TCTATGA		Eurogentec
TP53i9 minigene γ	CACCCTTCAGATGCTACTTG		Eurogentec
Neomycin	AGACAATCGGCTGCTCTGAT	CAATAGCAGCCAGTCCCT TC	Eurogentec
CPSF4	QuantiTect Primer Assay (product no 249900)		Qiagen
GAPDH	TGCACCACCAACTGCTTAG	G TTCAGCTCAGGGATGAC C	Eurogentec
TaqMan assays			
P53 (Exon8/9)	GAAGAGAATCTCCGCAAGA AAGG	TCCATCCAGTGGTTTCTT CTTG	Eurogentec
Probe P53 (Exon8/9)	AGCACTAAGCGAGCACTGCCCAACA		Eurogentec
p53 β	AACCACTGGATGGAGAATA TTTCAC	AAGAGAGCATGAAAATG GTTCTATGA	Eurogentec
Probe p53 β	CAGGACCAGACCAGCTTTCAA AAAAGAAAATTGTT		Eurogentec
Actin	GGCACCCAGCACAATGAAG	GCCGATCCACACGGAGT ACT	Eurogentec
Probe Actin	TCAAGATCATTGCTCCTCCTGAGCGC		Eurogentec