

SUPPLEMENTARY TABLE S2. FLOW CYTOMETRY PANEL FOR TIGIT STAINING

<i>Antibody</i>	<i>Fluorochrome</i>	<i>Clone</i>	<i>Amount used (μL/sample)</i>	<i>Catalog no.</i>	<i>Company</i>
CD95	BV421	DX2	5	562616	BD Biosciences
CD3	BV711	SP34-2	2.5	740807	BD Biosciences
CCR7	BV786	3D12	5	563710	BD Biosciences
CD4	PerCP-Cy TM 5.5	L200	10	552838	BD Biosciences
CD28	PE-CF594	CD28.2	5	562296	BD Biosciences
CD45	Alexa Fluor 700	D058-1283	5	561288	BD Biosciences
CD8a	Pe-Cy7	RPA-T8	2.5	557760	BD Biosciences
TIGIT	PerCP-eFluor 710	MBSA43	5	46-9500-42	eBioscience
LIVE/DEAD	Aqua	—	1.5	L34957	Invitrogen