

SUPPLEMENTARY TABLE S2. PRECISE II WITH UPDATED GLUCOSE CALCULATION ALGORITHM
AGREEMENT TO REFERENCE

<i>Reference glucose range (mg/dL)</i>	<i>No. of paired CGM system-YSI reference readings</i>	<i>Percent of CGM system readings within</i>			
		<i>15/15% of reference</i>	<i>20/20% of reference</i>	<i>30/30% of reference</i>	<i>40/40% of reference</i>
≤54	224	89	94	99	100
>54 to <71	848	87	95	98	99
71 to 180	8122	85	93	98	99
>180	8122	85	93	98	99
Total	15,753	87	94	99	100