

<b>A</b>	MS marker	Forward Primer	Reverse Primer	Tm	CEQ Fragment Size for NHP4026	CEQ Fragment Size for NF54	CEQ Fragment Size for NHP4076	CEQ Fragment Size for NHP4333	CEQ Fragment Size for NHP1337	CEQ Fragment Size for NHP4302	CEQ Fragment Size for NHP4373	CEQ Fragment Size for NHP3032
	<b>TA119</b>	TCCTCGATTATATTA TTGCA	TAATACATTCCCATT AGATG	53	223	241	235	229	230	238	230	230
	<b>TA81</b>	TGGACAAATGGGAA AGGAT	TTTCACACAACACA GGATT	53	135	120	126	129	122	119	116	116
	<b>TA77</b>	AAGCAAAAGAAAAT AAGGAA	GCACATAGTTGGTT CTCAC	53	197	197	197	-	-	197	197	206
	<b>TA62</b>	ATTATAGTGAAGAT AGCCA	TTACTTTCATCACTA TTTG	48	164	173	164	-	164	164	164	164

**B**

**NHP4026**

<b>NF54</b>	TA119	<b>NF54</b>										
<b>NHP4076</b>	TA119	TA119	<b>NHP4076</b>									
<b>NHP4333</b>	TA119	TA119	TA119	<b>NHP4333</b>								
<b>NHP1337</b>	TA119	TA119	TA119	TA81	<b>NHP1337</b>							
<b>NHP4302</b>	TA119	TA62	TA81	TA119	TA119	<b>NHP4302</b>						
<b>NHP4373</b>	TA119	TA119	TA81	TA81	TA81	TA119	<b>NHP4373</b>					
<b>NHP3032</b>	TA119	TA119	TA81	TA81	TA81	TA119	TA77					